

Restriction Digest and Electrophoresis without DNA Recovery

Background:

Restriction enzymes such as *EcoRI* recognize specific sequences of DNA and cleave the phosphodiester bond on each strand at that sequence. After digestion with a restriction enzyme, the resulting DNA fragments can be separated by agarose gel electrophoresis, and their corresponding sizes can be estimated.

In this experiment, we will perform a restriction digest using the restriction enzyme *EcoRI* on λ -phage DNA that has already been digested by the enzyme *HindIII*. (Please see the attached sequence site sheet for exact fragment lengths)

The recognition site of *EcoRI* is:

5'...G'AATTC...3'

3'...C'TTAAG...5'

Materials:

λ -DNA/*HindIII* solution (5 μ l = 2.5 μ g)

EcoRI enzyme (0.5 μ l = 10 units)

10X *EcoRI* buffer (5 μ l)

DD-H₂O (double-distilled)

37°C incubator

65°C bath

Micropipettes and tips

Gel dish

6X loading dye

1X TAE Buffer

1X TAE-agarose gel

100 bp dyed DNA ladder

Procedure:

In a one ml centrifuge tube add the following components in this order:

5 μ l λ -phase DNA/*HindIII*

5 μ l 10X *EcoRI* buffer

39.5 μ l DD-H₂O

0.5 μ l *EcoRI* enzyme


Mix the final solution with the pipette tip. *Do not vortex as this will rip the DNA.* Spin down the tube in the micro-centrifuge to ensure all liquid is at the bottom.

Place the tube in the 37° C bath for about 2 hours to facilitate digestion.

Meanwhile, clean the gel dish with ethanol and DI water. Heat the agarose gel to about 65 °C to liquefy, and pour about 35ml (1/3 full) into the sealed gel tray ensuring there are *no bubbles* in the liquid (ask for help here). Make sure the lane comb is properly in place before pouring the gel. Allow about 25 minutes for the gel to cool and solidify. Once fully solidified, remove the lane comb, and rotate the gel tray 90 degrees so that the gel wells are exposed to the negative terminal reservoir. Fill each of the reservoirs with 1X TAE buffer until the gel is submerged with about 1mm of buffer. Set the power supply to 60VDC and 120 minutes, but do not start it yet.

About now your digestion should be complete. Place the centrifuge tube into the 65°C bath for about 5 minutes. This effectively destroys the enzymes and stops digestion. Remove the centrifuge tube from the bath and let it cool for a few minutes. Place 10 μ l of 6X loading dye in the centrifuge tube and stir gently with the pipette tip. You are now ready to start the gel.

Pick someone in your group with a steady hand to place about 15 μl of your digest in the first lane of the gel, 15 μl of the dyed 100bp ladder solution in the third lane, and 15 μl of the dyed λ -phage DNA/*HindIII* solution in the fifth lane. Secure the electrophoretic lid, and start the power supply. Be careful not to shake or jostle the gel. Within a few minutes you should be able to see the loading dye moving in the gel, as shown below.



When the 120 minutes is up, the power supply will automatically turn off. Carefully remove the gel tray and gel from the TAE buffer solution and place it in the Ethidium Bromide staining bath for approximately 8 minutes. PLEASE BE CAREFUL: ETHIDIUM BROMIDE IS A DANGEROUS MUTAGEN AND CARCINOGEN, DO NOT SPLASH OR SPILL IT!! YOU MUST WEAR NITRILE GLOVES, SAFETY GLASSES AND A LAB COAT DURING THIS STEP. Make sure to thoroughly rinse the gel tray and spatula after contact with the Ethidium Bromide. Next rinse the gel in the water bath for 20 minutes. Finally, carefully place the gel on the UV light stage. Place the UV blocking cover over the stage and turn the UV lamp on. Quickly position the camera over the gel and take the best possible image, preferably one with a longer exposure. Admire your gel, you are now done.

Using a computer, you should be able to compare the calibration for the 100 bp ladder (found below) with your digested DNA, and verify the combined *HindIII* and *EcoRI* cut lengths. Does a linear relationship exist between the fragment length and the distance it goes in the gel?




Fig: 0.5 μg of 100 bp DNA Ladder visualized by ethidium bromide staining on a 1.3% TAE agarose gel (courtesy NEB)




Fig: Top lane- λ -phage DNA digested by *HindIII* & *EcoRI*, bottom lane: 100bp ladder. This is a decent experimental result from this protocol.

Lambda DNA: Location of Sites

Enzyme	#	Locations				
Apa I	1	10086				
Bmt I	1	34679				
BsiW I	1	19323				
Kas I	1	45679				
Nae I	1	20040				
Nar I	1	45679				
NgoM IV	1	20040				
Nhe I	1	34679				
PaeR7 I	1	33498				
PspOM I	1	10086				
SanD I (x)	1	28797				
Sfo I	1	45679				
SnaB I	1	12188				
Tii I	1	33498				
Xba I	1	24508				
Xho I	1	33498				
Acc65 I	2	17053	18556			
Afe I	2	20995	37057			
Asc I	2	3521	16648			
Avr I	2	24322	24396			
Bsa I	2	11424	42715			
Bsu36 I	2	26717	34318			
Eag I	2	19944	36654			
FspA I (x)	2	21804	21825			
Kpn I	2	17053	18556			
Pci I	2	628	39395			
PfiF I	2	11202	36120			
Pme I	2	8459	16293			
Sac I	2	24772	25877			
Sal I	2	32745	33244			
Tth111 I	2	11202	36120			
Afl II	3	6540	12618	42630		
BstZ17 I	3	15260	18834	19473		
Drd I	3	5116	9104	11090		
EcoO109 I	3	2815	28797	48473		
Pml I	3	26529	41482	42362		
PpuM I	3	2815	28797	48473		
Pvu I	3	11933	26254	35787		
Sma I	3	19397	31617	39888		
Xma I	3	19397	31617	39888		
ApaL I	4	5619	21798	27173	40216	
Bmr I	4	7054	11608	25691	30332	
Nco I	4	19329	23901	27868	44248	
SacI I	4	20320	20530	21606	40386	
BamH I	5	5505	22346	27972	34499	41732
BsrG I	5	5220	6142	15855	29392	32496
EcoR I	5	21226	26104	31747	39168	44972
Nru I	5	4590	28050	31703	32407	41808
RsrI I	5	3800	6041	13983	19288	22242
Sbf I	5	2555	2819	11834	19832	37000
Sca I	5	16421	18684	25685	27263	32802
SexA I	5	22264	31009	32838	40497	44408
Bgl II	6	415	22425	35711	38103	38754
Blp I	6	10297	10682	11661	16518	20744
BssH II	6	3522	4126	5627	14815	16649
SgrA I	6	7064	8680	12878	15653	16974

Enzyme	#	Locations				
Sph I	6	2212	12002	23942	24371	27374
Stu I	6	12434	31478	32997	39992	40596
Acl I	7	13529	16290	22580	22595	24642
Ban II	7	581	10086	19763	21570	24772
BbvC I	7	8012	18147	18465	30916	31222
BstB I	7	18048	25884	27980	29150	30396
Hind III	7	23130	25157	27479	36895	37459
Mlu I	7	458	5548	15372	17791	19996
Nde I	7	27630	29883	33679	36112	36668
PshA I	7	8920	9394	13512	15412	36925
Ava I	8	4720	19397	20999	27887	31617
Bcl I	8	8844	9361	13820	32729	37352
BsoB I	8	4720	19397	20999	27887	31617
BspH I	8	889	4650	4989	10249	18275
BssS I	8	20356	25572	27956	29425	34430
Mfe I	8	22687	22715	23054	25863	35764
Acc I	9	2190	15260	18834	19473	31301
Ahd I	9	6398	11238	12477	12915	16588
EcoN I	9	13509	21292	22377	25174	25223
Aat II	10	5105	9394	11243	14974	29036
Bae I	10	694	7665	13267	15775	16271
Bme1580 I	10	5619	5664	10086	11414	13039
Dra III	10	2954	5613	6635	8999	14477
Sap I	10	2397	6489	8702	10370	13286
Sty I	10	19329	21211	23901	24322	24396
Zra I	10	5105	9394	11243	14974	29036
Aar I (x)	12	554	10835	13941	14383	16164
Psi I	12	2285	9011	18941	19573	22014
Xcm I	12	958	4770	5065	5874	9193

There are no restriction sites for the following enzymes:

AsiS I, Fse I, I-Ceu I, I-Sce I, Not I, PI-Psp I, PI-Sce I, Pac I, Sfi I, Spe I, Srf I (x), Swa I

(x) = enzyme not available from NEB