

## Aph 162- Molecular Biology

### Day 3

#### Qiagen PCR purification of Double Digested YFP PCR product

Perform as on Day 1.

#### Nanodrop reading of purified, double digested YFP PCR product

See TA demo, as on Day 1.

#### Ligation of Vector and Insert

This step will incorporate our gene of interest into the open, linear vector we created during the double digest.

Reagent:	Reaction:	Control:
Double digested backbone DNA (from pZE21_LacZ)*	50 ng	50 ng
Insert DNA	Calc.	0 $\mu$ l
Quick Ligase Reaction Buffer (2X)	10 $\mu$ l	10 $\mu$ l
Quick Ligase	1 $\mu$ l	1 $\mu$ l
DDH <sub>2</sub> O	Calc.	Calc.
<b>Total Volume</b>	<b>21 <math>\mu</math>l</b>	<b>21 <math>\mu</math>l</b>

\* Your TA will provide the concentration.

#### Procedure:

- To maximize production of vector + insert ligated product, we add the insert in molar excess.
  - Calculate the amount (mass and corresponding volume) of insert DNA required to have a 1:3 molar ratio of vector to insert.**
    - Vector length - ~2200 bp
    - Insert length - ~720 bp
- To a microcentrifuge tube, add the appropriate amount of DNA and record the total volume. If the volume is below 10 $\mu$ l, add DDH<sub>2</sub>O to bring the volume to 10  $\mu$ l.
- Add 10  $\mu$ l of Quick Ligase Reaction Buffer and mix well.
- Ask your TA to add 1  $\mu$ l of the Quick Ligase, which contains an optimized ligation enzyme.
- Incubate at room temperature for 5 mins. Then put the reaction on ice, until **Cell Transformation** (see below).
- Perform a no-insert control by repeating these steps and omitting the insert – in its place add an appropriate amount of DDH<sub>2</sub>O.

**What does this control test for?**

#### PCR Purification of Ligated Product

Perform as on Day 1. Purify both the ligated product and the control.

#### Cell Transformation

A culture of the *E. coli* strain DH5 $\alpha$  (plasmid-free) was previously made competent, that is, able to accept a plasmid, following a protocol for electrocompetent cells.

*Note:* Make sure your ligated DNA has been purified before this step. See the included “PCR Purification Protocol” from Qiagen.

Procedure:

1. Warm up 1 mL LB media per transformation in a culture tube at 37°C.
2. Put two electroporation cuvettes and two vials of electrocompetent cells on ice.
  - a. Allow the tube of cells to thaw on ice.
1. Pipette 2  $\mu$ l (~8ng) of ligated DNA into the tube containing the cells and mix gently.
2. Pipette 50  $\mu$ l of this mixture into the electroporation cuvette, and electroporate at 1.8 kV. The time constant ( $\tau$  value) should be 3-4 msec.
3. Immediately pipette 750  $\mu$ l of warm LB from the culture tube into the cuvette, mix gently, and put back into the culture tube.
4. Incubate culture tube with cap loose on shaker at ~225 rpm for 1 hour at 37°C.
5. Plate culture on LB + Kanamycin + X-Gal plates
6. Let the plate dry, then incubate it upside-down overnight at 37°C.
7. With two separate transformations, make sure to transform cells with both the ligated product and the no-insert control.