

Bi 1x, Spring 2010


Week 6

- **Session 1**
 - TOPO cloning reaction
 - Bacterial transformation and plating

Background

TOPO mediated DNA cloning

When sequences are amplified with PCR, the Taq polymerase enzyme adds an extra (untemplated) deoxyadenosine ("A") on the 3'-end of the DNA fragments. These "A-overhangs" are complementary to the deoxythymidine "T" overhangs on the cloning site of the vector. The vector is supplied in the linear form and when you add a PCR product the complementary ends base pair to form a circular plasmid containing both vector and insert. The TOPO name in the kit refers to the fact that the vector contains topoisomerase I that is covalently bound to the vector at the insertion site. The topoisomerase enzyme has a ligase activity that ligates the phosphodiester backbone of the insert to the vector.¹



The vector

In addition to 5'- "T-overhangs", the plasmid vector pCR4-TOPO contains the genes encoding for ampicillin and kanamycin resistance. Bacteria transformed with the plasmid can then be selected with media that contains either antibiotic (in this case, we will use kanamycin because it exhibits greater stability). The pCR4-TOPO vector contains a lethal gene, *ccdB*, which is fused

to the LacZa gene. Upon disruption with a ligated PCR product, it permits only the growth of positive recombinants--i.e. cells that take up a vector without an insert will destroy themselves.

The cells

We will be using *E.coli* TOP10 competent cells ("One-Shot Cells") supplied with the kit. The term "competent" refers to the fact that these cells have been chemically treated to allow rapid uptake of DNA such as the pCR4-TOPO vector. The cells will be grown on agar plates containing Luria-Bertani broth (LB) and kanamycin.

Purification:

Perform PCR Purification from Qiagen (see attached)

Elute 40 ul in water

Preparation of overhangs:

Prepare on ice in a PCR tube:

8uL purified PCR

1uL 10X PCR buffer (from TOPO kit)

0.5uL dNTP (from TOPO kit)

0.5uL Taq (from TOPO kit)

10 μ l total volume

Incubate at 72° C for 20min. in PCR machine

Cloning reaction

1. Set up the following reaction on ice:

4 μ l	PCR product
1 μ l	salt solution
1 μ l	TOPO vector*

6 μ l	Total volume

*TOPO vector is extremely expensive. A TA will add the TOPO vector to your reaction when you are ready.

2. Mix gently with a pipettor and incubate at room temperature for 10 minutes. After the reaction is complete, place the tube back on ice.

During this brief incubation time, the topoisomerase enzyme will ligate the ends of the pCR4-TOPO vector to your PCR products.

Transformation reaction & plating

1. Place a tube of the TOP-10 competent cells on ice to thaw. Five minutes after the tube is thawed, proceed to step 2.
2. Add 2 µl of the TOPO cloning reaction and mix gently into the cells by flicking the tube lightly with your finger several times. **Do not mix by pipetting up and down, as the cells are very delicate.**
3. Incubate the cells on ice for 30 minutes.
4. **Heat shock** the cells by placing the tube in a 42°C water bath for **exactly** 30 seconds.

Heat shocking causes these cells to uptake DNA by a mechanism that is not entirely understood. One theory is that the heat shock lowers the membrane potential to allow DNA to enter the cell more easily.²

5. Immediately transfer the tube to ice and add 250 µl of room temperature SOC medium.

SOC is a very rich medium that will help the cells recover after their “shocking” experience.

6. Ensure that the tube is screwed on tightly. Tape the tube horizontally to the bottom of the 37°C incubator. Incubate at 37°C, 200 rpm shaking for 30 minutes to allow the cells to recover.
7. Spread 50 µl of cells on one plate and 100 µl of cells on another plate. *See the demo performed by your TAs.* Incubate overnight at 37°C.

² Panja S, *et al.* “**How does plasmid DNA penetrate cell membranes in artificial transformation process of Escherichia coli?**” *Mol Membr Biol.* 2008 Aug;25(5):411-22.