

Bi 1X, Spring 2012

Week 2

- Session 1: Growth Curves
-

Background:

Cell challenge with anti-bacterial substances and Spectroscopy


While most bacteria pose no threat to humans, there are a number of species that can cause very serious infections. One major challenge facing biomedical scientists is how to kill off these bacteria after an infection has taken hold, a challenge that has become more and more challenging as the abuse of antibiotics has led to the widespread occurrence of antibiotic resistant strains. To better understand the concepts involved in combating bacterial infection, we will measure the bulk growth rate of *E. coli* in media after challenge with various anti-bacterial substances. We measure cell growth using spectroscopy and optical density as outlined in the following section. The included formulas will prove useful in analyzing your spectroscopic data.

Beer's Law and Optical Density

It is often extremely useful to know the concentration of cells in a liquid culture, and can be used in various kinetic measurements of fitness and protein production. Light scattering is the primary method by which we measure cell density. In the following section, we will discuss the basic physical concepts that allow determination of concentration as well as the formulas that relate optical properties to concentration.

The basic physics are as follows:

1. Light, with a wavelength comparable to the size of the object being measured (600 nm), enters the sample at an initial intensity (I_0)
2. Cells, having a different index of refraction as the surrounding medium, randomly reflect and scatter light out of the incident light path (scattering with no change in momentum). The amount of scattering is proportional to the amount of cells at any given point in the sample. (Caveat: to get an accurate reading, cell density must be low enough that light is only scattered *once* in its journey through the sample. Often, this means that cells must be serially diluted to get an accurate reading.)
3. Upon exiting the sample, the light intensity is reduced by $I_1 = I_0 e^{-\alpha lc}$, where α is the wavelength-dependent molar absorption coefficient of the species in question, l is path length – usually 1 cm in modern spectrophotometers, and c is the concentration, this is known as *Beer's Law*. See the figure below:



4. The transmission efficiency of the sample is the ratio of the final and initial light intensity, $T = I_t/I_0$, where the measured quantity, called “absorbance” is $A = \ln(I_0/I_t) = -\ln(T)$. The general rule of thumb is that accurate absorbance readings lie in the range $0.01 < A < 1$. (Serial dilutions of the sample will be necessary to bring A within this range.)

Your Mission:

- Measure baseline growth rate in rich media.
- Measure how various anti-bacterial substances influence actively growing E Coli.
- Determine the correspondence between OD₆₀₀ and cell density.

Protocol:

Prior to your involvement:

- 1) Groups of two will be assigned a particular mixture in which to assess growth rates.
 - i) Growth 1-6: NZC media
- 2) Overnight cultures of wild type (C600) *E. coli* were inoculated into 50mL of each of these media, and incubated on a shaker at 37°C for 15 minutes prior to the start of the experiment. This ensures that the so-called “lag phase” has passed, and that (slow) exponential growth has begun.

Procedure:

1. Every 10-15 minutes remove a 1mL sample of your culture using a sterile pipette. Deposit this into a 1.5mL cuvette and immediately bring to the spectrophotometer (the “spec”).
2. Carefully insert the cuvette into the spec in the correct orientation. Make sure the spec is set to read OD₆₀₀, and measure the absorbance. Record the reading and the EXACT clock time.

NOTE: You should have been given a “blank”, i.e. a sealed cuvette with sterile media. Absorbance is always measured relative to the blank. **Be sure to measure the absorbance of the blank before every culture measurement!**

3. Properly dispose of the sample cuvette, and repeat these steps for the next 45 minutes.
4. After the first hour is up, take one more measurement, and then add your assigned amendment.
 - i. Growth 2: Add nothing
 - ii. Growth 2: Add 3×10^{-6} moles Ampicillin (263uL of stock)
 - iii. Growth 3: Add 3×10^{-6} moles Kanamycin (148uL of stock)
 - iv. Growth 4: Add 3×10^{-6} moles Rifampicin (250uL of stock)
 - v. Growth 5: Add 10uL T7 stock (10^6 pfu/uL)
 - vi. Growth 6: Add 10uL T7 stock (10^6 pfu/uL) and 3×10^{-6} moles Kanamycin
5. Take the last cuvette that you OD'd and follow the next few steps.

